



Isolation of Bacteria Using the Total RNA Kit I (R6834)

- 1. Grow Bacteria in LB media to log phase. (Do not use an overnight culture.)**
2. Harvest no more than 3 ml culture ($<1 \times 10^9$ bacteria) by centrifugation at 4,000-5000 x g for 5- min at 4°C.
3. Resuspend the Bacterial Pellet with the appropriate amount of TE Buffer Containing Lysozyme (1 mg/mL for Gram Negative, 4 mg/mL for Gram Positive). Incubate at room temperature for 4 minutes.
3. Discard medium and resuspend cells in the appropriate amount of TRK Lysis Buffer according the chart below. Mix by pipetting up and down 5-10 times.

Number of Cells	Amt of TE Buffer(containing Lysozyme)	Amt of EtOH (96-100%)	Amount of TRK Lysis Buffer
$<5 \times 10^8$	100 μ l	250 μ l	350 μ l
$5 \times 10^8 - 1 \times 10^9$	200 μ l	500 μ l	700 μ l

4. Add the appropriate amount of EtOH to the sample according to the chart above. Mix by Vortexing.
5. Follow the Total RNA Kit I protocol from Step 6 of the Animal Cells Protocol