

Protocol for Isolating RNA from Human Blood Collected into PAXgene™ Blood RNA Tube

The following modified protocol is designed for isolating RNA from human blood samples collected with PAXgene™ RNA tube using E.Z.N.A. Blood RNA Kit (R6814-01/02).

1. Centrifuge the PAXgene™ Blood RNA Tube for 10 minutes at 3000-5000 x g using a swing-bucket rotor.
2. Remove the supernatant by decanting or pipetting. Add 4 ml RNase-free water to the pellet, close the tube with a cap.
3. Vortex until the pellet is completely resuspended. Centrifuge at 3000-5000 x g for 10 minutes using a swing-bucket rotor. Remove and discard the entire supernatant. **Note: Incompletely removal of supernatant will reduce the lysis efficiency and dilute the lysate, and therefore reduce the RNA yield.**
4. Add 300µl TRK Lysis Buffer, vortex the sample until the pellet is completely dissolved.
5. Transfer the sample into a new 1.5 ml microcentrifuge tube, add 590µl of DEPC-water and 10 µl Proteinase K (25mg/ml). Mix thoroughly by vortexing for 5 seconds. Incubate at 55°C for 10 minutes using a shaker-incubator.
6. Transfer the sample directly into a Homogenizer column placed into a 2 ml collection tube, centrifuge at maximum speed (>13,000 x g) for 3 minutes.
7. Carefully transfer the entire supernatant of the flow-through fraction to a new 1.5 ml tube without disturbing the pellet in the 2ml collection tube.
8. Add 450µl of absolute ethanol (96-100%). Mix the sample thoroughly by vortexing.
9. Pipet 750µl of sample into the HiBind RNA column placed in a 2 ml collection tube. Centrifuge at 8000-20,000 x g for 1 minute. Discard the flow-through and collection tube.
10. Place the HiBind RNA column into a new 2 ml collection tube. Transfer the remaining sample into the HiBind RNA column with pipettor. Centrifuge at 8000-20,000 x g for 1 minute. Discard the flow-through and re-use the collection tube.

Note: this is the starting point for optional On-Membrane DNase I digestion procedure. See detail information from standard user manual. If the optional DNase I digestion is desired, using step 11a, otherwise, go to step 11b.

- 11a. Place the HiBind RNA column into 2 ml collection tube and add 350 μ RNA wash Buffer I into the HiBind RNA column. Spin at 8000-20,000 x g for 1 minute. Continue with Optional DNase I digestion protocol. (See standard user manual for details)
- 11b. Add 500 μ RNA wash Buffer I into the HiBind RNA column and spin at 8000-20,000 x g for 1 minute. Discard the flow-through and re-use the collection tube.
12. Place the HiBind RNA column into same collection tube. Add 600 μ RNA wash Buffer II into the HiBind RNA column and spin at 8000-20,000 x g for 1 minute. Discard the flow-through and re-use the collection tube.
13. Place the HiBind RNA column into same collection tube. Wash the column again by adding another 600 μ RNA wash Buffer II into the HiBind RNA column and spin at 8000-20,000 x g for 1 minute. Discard the flow-through and re-use the collection tube.
14. Place the HiBind RNA column into the same collection tube and centrifuge at maximum speed for 2 minutes.
15. Add 50-70 μ l DEPC-water or RNase-Free water directly onto the center of the membrane in the HiBind RNA column. Incubate 1 minute at room temperature. Centrifuge at maximum speed for 2 minutes to elute RNA.

Paxgene™ is a trademark of Qiagen N.V.