

DNA Ladders

Omega Biotek offers a wide range of products for the accurate size and mass determination of DNA. Most of DNA Ladders includes one or two higher-intensity bands for easy identification and orientation. The ladders are supplied premixed with loading buffer and are stable at room temperature for a few weeks. There is no need to heat, mix, or dilute ladders prior to loading them on gel.

50bp DNA Ladder plus

- CAT.No.: M07-01
- Packing: 250ul
- Concentration: 0.128 ug/ul
- Store: at 4°C for periods up to 6 months. For longer periods store at -20°C.
- Description: 50 DNA Ladder plus is designed for sizing and approximate quantification of wide range linear double-stranded DNA fragments on agarose gel. The ladder is composed of thirty purified individual DNA fragments (in base pairs): 1000, 900, 800, 700, 600, 500, 400, 300, 250, 200, 150, 100, 50. It contains two high-intensity bands (500bp and 250 bp) for easy identification and orientation.
- Direction:
 - ◇ This product can load directly, since it contains 1× loading buffer. It is free of fluorescence dye such as EtBr.
 - ◇ For a loading well less than 6mm in width, 5ul of this product is enough for clear bands. For a wider well, more product may be required.
 - ◇ Gel and buffer may effect the electrophoresis. For better result, please use high quality agarose and fresh running buffer. 3.0% agarose gel is recommended.
- Electrophoresis Map:

