

Contents

| | |
|---|----|
| Introduction | 2 |
| Kit Contents | 3 |
| Before Starting | 3 |
| Protocol 1: Fastfilter Endo-Free Plasmid Midi Kit Spin Protocol | 4 |
| ■ Growth of bacterial culture | 4 |
| ■ Lyse bacterial cells with Alkaline-SDS Solution | 4 |
| ■ Clear the lysate with Lysate Clearance Filter Syringe | 5 |
| ■ Remove Endotoxins with ETR Binding Buffer | 5 |
| ■ Purify Plasmid DNA with HiBind™ DNA Midi Column | 6 |
| ■ Elution Plasmid DNA From HiBind™ DNA Midi Column | 7 |
| ■ Alternative Protocol of Elution Plasmid | 7 |
| Protocol 2: Fastfilter Endo-Free Plasmid Midi Kit Vacuum Protocol | 9 |
| Protocol 3: Low Copy-Number Plasmids Protocol | 9 |
| Yield and Quality of DNA | 10 |
| Plasmid Copy Number and Expected Yield | 10 |
| Trouble Shooting Guide | 11 |
| Ordering information | 12 |

Revised March 2009

Introduction

The E.Z.N.A.® family of products is an innovative system that radically simplifies extraction and purification of nucleic acids from a variety of sources. Key to the system is the Omega Bio-Tek's (OBI) proprietary HiBind™ matrix that avidly, but reversibly, binds DNA or RNA under certain optimal conditions allowing proteins and other contaminants to be removed. Nucleic acids are easily eluted with deionized water or low salt buffer.

Endotoxins are lipopolysaccharides (LPS), found in the outer cell membrane of the gram negative bacteria such as *E. coli*. In Mammalia system, the endotoxins are pyrogenic, it can cause fever and endotoxin shock syndrome. The endotoxin contamination is one of the major concern for gene therapy. The endotoxin contamination in plasmid can cause lower transfection efficiency for endotoxin sensitive cell lines. Since Endotoxin is negative charged molecular like DNA, both DNA and endotoxin moleculars behave similarly on the surface of silica and anion-exchange chromatography which are the most common used technologies for plasmid purification.

The E.Z.N.A.® Fastfilter Endo-Free Plasmid Midi Kit combines the power of HiBind™ technology with the time-tested consistency of alkaline-SDS lysis of bacterial cells to deliver high quality DNA. Omega Bio-Tek's HiBind™ Midi columns facilitate the binding, washing, and elution steps thus enabling multiple samples to be simultaneously processed. This kit also include a special filter cartridge, which replaces the centrifugation step following alkaline lysis. Following lysis of the cell, endotoxins are removed from the cleared cell lysate with simple extraction-phase-separation steps. Then DNA is bound to the silica membrane and contaminants are removed with a simple wash step. Yields vary according to plasmid copy number, *E. coli* strain, and conditions of growth, but up to 100-250 µg of high copy-number plasmid or 10-100 µg of low copy number plasmid can be purified from 30-50 ml overnight culture. Up to 100 ml bacterial cultures can be used when working with low copy number vector. The product is suitable for automated fluorescent DNA sequencing, restriction endonuclease digestion, transfection of mammalian cells, and other manipulations.

Storage and Stability

All E.Z.N.A.® Plasmid isolation components are guaranteed for at least 12 months from the date of purchase when stored as follows: Solution I (once RNase A is added) and ETR Binding Buffer at 4°C, all other material at 22-25°C.

Kit Contents

| Product No. | D6915-00B | D6915-01B | D6915-03B | D6915-04B |
|---------------------------------|-----------|-----------|-----------|------------|
| Purification Times | 2 preps | 10 Preps | 25 Preps | 100 Preps |
| HiBind™ DNA Midi Columns | 2 | 10 | 25 | 100 |
| Lysate Clearance Filter Syringe | 2 | 10 | 25 | 100 |
| Solution I | 6 ml | 30 ml | 80 ml | 270 ml |
| Solution II | 6 ml | 30 ml | 80 ml | 270 ml |
| Buffer N3 | 4 ml | 20 ml | 40 ml | 140 ml |
| ETR Binding Buffer | 15 ml | 60 ml | 170 ml | 650 ml |
| ETR Wash Buffer | 5 ml | 30 ml | 60 ml | 230 ml |
| Buffer EHB | 8 ml | 35 ml | 85 ml | 2 x 170 ml |
| DNA Wash Buffer | 4 ml | 20 ml | 50 ml | 3 x 50 ml |
| RNase A | pre-added | 150 µl | 400 µl | 1.4 ml |
| Endotoxin-Free Elution Buffer | 8 ml | 20 ml | 60 ml | 2 x 120 ml |
| Instruction Booklet | 1 | 1 | 1 | 1 |

Before Starting

Briefly examine this booklet and become familiar with each step. Prepare all components and have the necessary materials ready before starting.

Supplied By Laboratory centrifuge equipped with **swinging-bucket** rotor
User: capable of 3,000-5,000 x g.
 Sterile 15 ml & 50 ml centrifuge tubes capable of 5000 x g.
 96-100% ethanol

IMPORTANT

1. Add vial of **RNase A** to the bottle of **Solution I** provided and store at 4°C.
2. **DNA Wash Buffer** is to be diluted with **96-100% ethanol** as follows:

| | |
|-----------|---|
| D6915-00B | Add 16 ml of absolute ethanol to bottle |
| D6915-01B | Add 80 ml of absolute ethanol to bottle |
| D6915-03B | Add 200 ml of absolute ethanol to bottle |
| D6915-04B | Add 200 ml of absolute ethanol per bottle |

Store diluted DNA Wash Buffer at room temperature !

Protocol 1: Endo-Free Plasmid Midi Kit Spin Protocol

This Protocol is designed to isolate 100-250 µg of high Copy-Number plasmids or 10-100 µg of low Copy-Number Plasmid from 30-50 ml overnight cultures using E.Z.N.A.® Fastfilter Endo-Free Plasmid Midi Kit. For increasing yield of low Copy-Number plasmid, proceed as “Low Copy-Number Plasmids protocol” on page 8.

■ Growth of bacterial culture

1. **Inoculate 20-50 ml LB/ampicillin (50 µg/ml) medium placed in a 200-400 milliliter culture flask with *E.coli* carrying desired plasmid and grow at 37°C with agitation for 12-16 h.** For best results use overnight culture (0.3-0.5 ml) as the inoculum. It is strongly recommended that an *endA* negative strain of *E.coli* be used for routine plasmid isolation. Examples of such strains include DH5α® and JM109®.

Optimal growth conditions of bacteria is vital of obtaining maximal plasmid DNA yields. The best conditions are achieved by picking a single isolated colony from a freshly transformed or freshly plate to inoculate a 2-5ml starter culture containing the appropriate antibiotic. Incubate for ~8hr at 37°C with vigorous shaking (~300rpm). Then used to inoculate appropriate volume of Pre-warmed liquid growth medium with antibiotic. Grow at 37°C for 12-16 hr with vigorous shaking(~300rpm).Using a flask or vessel with a volume of a least 3-4 times the volume of the culture and dilute the starter culture 1/500 to 1/1000 into growth medium.

Following overnight bacterial growth, an OD₆₀₀ of 1.5-2.0 indicates a well-grown culture. For the best result determination of OD₆₀₀ for each culture is recommended. it is important to dilute the bacterial culture (10 to 20 fold) to enable photometric measurement in the linear range between 0.1 and 0.5 OD₆₀₀. We recommend a bacterial density of between 2.0 and 3.0 at OD₆₀₀. When using untrient-rich media, care should be taken ensure that the cell density does not exceed an OD₆₀₀ of 3.0.

If using a frozen glycerol stock as inoculum, streak it onto an agar plate containing the appropriate antibiotic for single colony isolation. Then picking a single colony and inoculate the 2-5ml starter culture as described above.

■ Lyse bacterial cells with Alkaline-SDS Solution

2. **Pellet up to 20-50 ml bacteria in appropriate vessels by centrifugation at 3,500-5,000 x g for 10 min at room temperature.**
3. Decant or aspirate medium and discard. To ensure that all traces of the medium are removed, use a clean paper towel to blot excess liquid from the wall of the vessel. **To the bacterial pellet add 2.5 ml Solution I/RNase A.**

Resuspend cells completely by vortexing or pipetting up and down. *Complete resuspension of cell pellet is vital for obtaining good yield.*

4. **Add 2.5 ml Solution II and mix gently but thoroughly by inverting and rotating tube 7-10 times to obtain a cleared lysate.** A 3-5 min incubation at room temperature may be necessary.

Avoid vigorous mixing as this will shear chromosomal DNA and lower plasmid purity. (Store Solution II tightly capped when not in use.)

5. Prepare a Lysate Clearance Filter Syringe by removing the plunger and place the barrel in a tube rack to keep the syringe upright.
6. **Add 1.25 ml Buffer N3 and gently mix by inverting tube several times until a flocculent white precipitate forms.** This may require a 2-3 min incubation at room temperature with occasional mixing. **Prepare a Lysate Clearance Filter Syringe by placing the barrel in a tube rack to keep the syringe upright.**

The Buffers must be mixed thoroughly. If the mixture appears still viscous, brownish and conglobated, more mixing is required to completely neutralize the solution. Complete neutralization of the solution is vital of obtaining good yields. Using ice-cold Neutralization Buffer is helpful to precipitate more bacterial proteins.

7. **Prepare the HiBind Midi Column.** Place a HiBind Midi column into a 15 ml collection tube, provided. Add 2 ml of Buffer GPS to the column and let it sit at room temperature for 3-10 min. Spin in a swinging bucket rotor at 3,000-5,000xg for 5 min at room temperature. Discard the eluate and assemble the column again in the 15 ml collection tube.

■ **Clear the lysate with Lysate Clearance Filter Syringe**

8. **Immediately pour the lysate into the barrel of the Lysate Clearance Filter Syringe. Allow the cell lysate to sit for 2 minutes.** The white precipitate should float to the top. Use a new 50 ml tube to collect the cell lysate. Insert the plunger back into the barrel of the syringe.
9. **Hold the Lysate Clearance filter syringe barrel over the 50 ml tube and gently insert the plunger to expel the cleared lysate to the tube.**

Note: Some of the lysate may remain in the flocculent precipitate, do not force this residual lysate through the filter. Do not use a new 15 ml tube to collect the cell cleared lysate, because ETR Solution may not be removed in step 11 with too much solution in 15 ml tube.

Alternatively, the cell debris and KDS-precipitation can be removed by

centrifugation at 15,000 x g for 10 min at 4°C, instead of using Clearance Filter Syringe in step 8-9. A tightly packed cell debris pellet indicates efficient lysis.

10. **Add Equal volume of ETR Binding Buffer to the filtered lysate, mix by inverting the tube 7-10 times.**
11. **Add 3.5 ml of the lysate from step 10 into the HiBind® DNA Midi column assemble into the 15 ml Collection tube, centrifuge at 3,000-5000 x g for 3-5 min.** Discard the flow-through and re-use the collection tube.
12. **Repeat step 11 until all the remaining of the cleared lysate pass through the HiBind® DNA Midi Column**
13. **Add 2.0 ml ETR Wash Buffer to the DNA Midi column and centrifuge as above.** Discard the flow-through and re-use the collection tube.
14. **Add 3.0 ml Buffer EHB to the DNA Midi column and centrifuge as above.** Discard the flow-through and re-use the collection tube.
15. **Add 3.5 ml DNA Wash Buffer diluted with absolute ethanol to the DNA Midi column and centrifuge as above.** Discard the flow-through and re-use the collection tube.

Note: DNA Wash Buffer Concentrate must be diluted with absolute ethanol before use. See label for directions. If refrigerated, DNA Wash Buffer must be brought to room temperature before use.

16. **Add 3.5 ml DNA Wash Buffer diluted with absolute ethanol to the DNA Midi column and centrifuge as above.** Discard the flow-through and re-use the collection tube.
 17. **Centrifuge the empty column at maxi speed (no more than 8,000 x g) for 10-15 min to dry the column matrix.** Do not skip this step - it is critical for removing ethanol from the column.
- **Elution Plasmid DNA From HiBind™ DNA Midi column.**
Optional: For maximal yield and high concentration of plasmid, see alternative protocol of elution. For fast elution, proceed step 18-19.
18. **Drying the column (Optional):** choose either of the methods below to further dry the column before eluting DNA.

1. Place the column into a vacuum chamber at room temperature. Any device connected to a vacuum source may be used. Seal the chamber and apply

vacuum for 15 min. Remove the column and proceed to step 20.

2. Bake the midi column at 65°C in a vacuum oven or incubator for 10-15 min. Proceed to step 20.
19. Place column into a clean 15 ml centrifuge tube. Add 0.5-1 ml (depending on desired final concentration) Endotoxin-Free Elution Buffer onto the column matrix and let it sit at room temperature for 2-3 min. Centrifuge at maxi speed (no more than 8000 x g) for 5 min to elute DNA. This represents approximately 70-80% of bound DNA. An optional second elution will yield any residual DNA, though at a lower concentration. Alternatively, a second elution may be performed using the first eluate to maintain a high DNA concentration. Preheating the water to 65°C and allowing the column to soak 2 min at room temperature before elution may significantly increase yields.

Note: The plasmid DNA obtained using this protocol performs well in PCR, restriction digests, lipid mediated transfection and transformation. The expected concentration of plasmid is vary between different copy number vector. However, the concentration of high copy-number plasmid is 150-400ug/ml. Some residual ethanol is present, but does not interfere with these downstream applications. One may get high concentration and absolutely remove ethanol with optional elution step as following.

■ Alternative Protocol of Elution Plasmid

1. Place HiBind™ DNA Midi column into a clean 15 ml centrifuge tube. Add 3 ml Endotoxin-Free Elution Buffer (or TE buffer) onto the column matrix and centrifuge at maxi speed (no more than 8,000 x g) for 2 min to elute DNA. Preheating the water to 70°C and allowing the column to soak 2 min at room temperature before elution may significantly increase yields.
2. Carefully transfer eluted plasmid from 15 ml centrifuge tube to a clean tube suitable for precipitation. Add 130 µl 5M NaCl and 2.2 ml room temperature isopropanol. Vortex to mix and centrifuge at >15,000 x g for 30 min at 4°C. Carefully decant the supernatant.
3. Wash DNA pellet once with 1ml ice-cold 70% ethanol and centrifuge at > 15,000 x g for 10 min. Carefully decant the supernatant without disturbing the pellet. Air-dry the pellet for 5-10 min.
4. Finally resuspend DNA pellet in 200-500 µl (depending on desired concentration of final product) Endotoxin-Free Elution Buffer or TE Buffer.

Protocol 2: Endo-Free Plasmid Midi Kit Vacuum Protocol

1. Prepare the cleared lysate by following step 1-12 of spin procedure on page 4-6.
2. Prepare the vacuum manifold according to manufacturer's instructions and connect the Midi column to the manifold.
3. Transfer the cleared lysate into the HiBind® DNA Midi column, be careful not to overfill the column, apply the vacuum to allow all sample pass through the column.
4. Add 2 ml ETR Wash Buffer to the column and apply the vacuum to draw the liquid through the column.
5. Add 3.5 ml EHB buffer to the column and apply the vacuum to draw the liquid through the column.
6. Wash the column: Add 3.5 ml of DNA Wash Buffer (pre-diluted with absolute ethanol) into the column and allow it pass through the column.
7. Wash the column again with 3.5 ml DNA Wash Buffer by repeating step 5. Keep the vacuum on for another 10-15 min after the liquid pass through the column. (This step to ensure the removal of residue ethanol).
8. Proceed step 18-19 of Spin protocol on page 7.

Protocol 3: Low Copy-Number Plasmids Protocol

Low copy number plasmids generally give 0.1-1 µg DNA per ml overnight culture. For isolation of plasmid from low copy-number plasmids (0.1-1 µg/ml culture) or low-midi copy-number plasmids (1-2 µg/ml culture) bacteria, the method can be modified to essentially increase the yield if necessary. Start with 50-100 ml bacterial culture, centrifuge for 10 min at 3,500-5,000 xg in a centrifuge tube. Proceed to Step 3 (Page 4) and double the volumes of Solution I, II, and Buffer N3. Continue as above using only one HiBind® DNA Midi column. There is no need to increase the volumes of Buffer HB and DNA Wash Buffer used.

Note: This method is not recommended for high copy number plasmids because above 50 ml culture, the HiBind® DNA Midi column quickly becomes saturated. In this situation we recommend processing of multiple samples from the same culture.

Yield and Quality of DNA

Determine the absorbance of an appropriate dilution (20- to 50-fold) of the sample at 260 nm and then at 280 nm. The DNA concentration is calculated as follows:

$$\text{DNA concentration} = \text{Absorbance}_{260} \times 50 \times (\text{Dilution Factor}) \mu\text{g/ml}$$

A value greater than 1.8 indicates greater than 90% nucleic acid. Alternatively, quantity (as well as quality) can sometimes best be determined by agarose gel/ethidium bromide electrophoresis by comparison to DNA samples of known concentrations. Typically, the majority of the DNA eluted is in monomeric supercoil form, though concatamers may also be present.

Plasmid Copy-Number and Expected Yield

The Yield and quality of the plasmid DNA depend on a number of factors including plasmid copy number, size of insert, host strain, culture volume, culture medium and binding capacity of the kit. In these factors, The copy number of vector, culture volume and binding capacity of the kit are most important. Copy number of plasmid is vary from one copy to several hundred copies per cell as dicated by their origin of replication. But very large plasmids are often maintained at very low copy number per cell. The expected yield of 200 ml overnight cultures (LB medium) are indicated in the table.

| Plasmid | Replicon | Copy Number | expected Yield of 50ml culture |
|----------------------------|----------|-------------|--------------------------------|
| pUC vector | pMB1 | 500-700 | 150-250 µg |
| pBR322 and its derivatives | pMB1 | 15-20 | 10-20 µg |
| PACYC and its derivatives | p15A | 10-12 | 5-10 µg |
| pSC101 and its derivatives | pSC101 | ~5 | 5 µg |
| pBluescript | ColE14 | 300-500 | 100-180 µg |
| ColE14 | ColE14 | 15-20 | 10-20 µg |
| pGEM | pMB1 | 300-700 | 100-200 µg |

Trouble Shooting Guide

| Problem | Likely Cause | Suggestions |
|--|---|---|
| Low DNA yields | Poor cell lysis | <p>Only use LB or YT medium containing ampicillin. Do not use more than 100 ml.</p> <p>Cells may not be dispersed adequately prior to addition of Solution II. Vortex/ pipet cell suspension to completely disperse bacterial clumps.</p> <p>Increase incubation time with Solution II to obtain a clear lysate.</p> <p>Solution II if not tightly closed, may need to be replaced. Prepare as follows: 0.2 N NaOH, 1% SDS.</p> |
| | Bacterial culture overgrown or not fresh. | Do not incubate cultures for more than 16 hr at 37°C. Storage of cultures for extended periods prior to plasmid isolation is detrimental to yield and quality. |
| | Low copy-number plasmid used | Such plasmids may yield as little as 0.1µg DNA from a 1 ml overnight culture. Increase culture volume to 100 ml. |
| No DNA eluted. | DNA Wash Buffer Concentrate not diluted with absolute ethanol. | Prepare Wash Buffer Concentrate as instructed on the label. |
| High molecular weight DNA contamination of product. | Over mixing of cell lysate upon addition of Solution II. | Do not vortex or mix aggressively after adding Solution II. Adequate mixing is obtained by simply inverting and rotating tube to cover walls with viscous lysate. |
| Optical densities do not agree with DNA yield on agarose gel. | Trace contaminants eluted from column increase A_{260} . | Make sure to wash column as instructed. Alternatively, rely on agarose gel/ethidium bromide electrophoresis for quantitation. |
| RNA visible on agarose gel. | RNase A not added to Solution I. | Add 1 vial of RNase to each bottle of Solution I. |
| Plasmid DNA floats out of well while loading agarose gel, does not freeze, or smells of ethanol. | Ethanol traces not completely removed from column following wash steps. | Centrifuge column at 3000 x g for 10 min to dry the column. A swinging-bucket rotor is recommended for centrifugation. Alternatively, precipitate the eluted DNA with isopropanol as indicated. |

Ordering Information

| Product | Applications | Cat. No. |
|--|--|----------|
| Standard E.Z.N.A.™ Plasmid Isolation System | | |
| Plasmid Mini Kit I | Isolation of up to 30 µg Plasmid DNA | D6942/3 |
| Plasmid Mini Kit II | Isolation of up to 70 µg Plasmid DNA | D6945 |
| Plasmid Midi Kit | Isolation of up to 250 µg Plasmid DNA | D6904 |
| Plasmid Maxip Kit | Isolation of up to 1.5 mg Plasmid DNA | D6922 |
| Fastfilter Plasmid Midi Kit | Isolation of up to 250 µg Plasmid DNA, featuring filter syringes for lysate clearance | D6905 |
| Fastfilter Plasmid Maxi Kit | Isolation of up to 1.5 mg Plasmid DNA, featuring filter syringes for lysate clearance | D6924 |
| E-Z 96 Fastfilter Plasmid Kit | Isolation of Plasmid DNA using a 96-well format | D1097 |
| E-Z 96 SE Plasmid Kit | Isolation of plasmid DNA using a single plate | D1095 |
| Yeast Plasmid Isolation Kit | Isolation fo Yeast Plasmid DNA | D3476 |
| E.Z.N.A.™ Endotoxin Free Plasmid Isolation System | | |
| Endo-Free Plasmid Mini Kit I | Isolation of up to 30 µg Endotoxin free Plasmid DNA | D6948 |
| Endo-Free Plasmid Mini Kit II | Isolation of up to 70 µg Endotoxin free Plasmid DNA | D6950 |
| Endo-Free Plasmid Mid Kit | Isolation of up to 250 µg Endotoxin free Plasmid DNA, featuring filter syringes for lysate clearance | D6915 |
| Endo-Free Plasmid Maxi Kit | Isolation of up to 1.5 mg Endotoxin free Plasmid DNA, featuring filter syringes for lysate clearance | D6926 |
| E.Z.N.A.™ H P Plasmid Isolation System | | |
| HP Plasmid Mini Kit I | Isolation of up to 30 µg of High Purity Plasmid DNA | D7042 |
| HP Plasmid MidiKit | Isolation of up to 200 µg of High Purity Plasmid DNA | D7004 |
| HP Plasmid Maxi Kit | Isolation of up to 1.5 mg of High Purity Plasmid DNA | D7022 |
| E.Z.N.A.™ Single Strand Phage DNA Isolation Kits | | |
| M13 Isolation Kit | Isolation of up to 15µg of single stranded phage DNA | D6900 |
| E-Z 96 M13 Isolation Kit | Isolation of up to 15µg of M-13 DNA using a 96-well format | D1900 |
| E.Z.N.A.™ Large Construct DNA Isolation Kits | | |
| BAC/PAC DNA Isolation Kit | Effective purification of BAC or PAC DNA | D2156 |
| BAC/PAC DNA Isolation Kit | Parallel purification of BAC or PAC using a 96-well format | D1056 |